

Product datasheet for **MR200624A1V**

Mouse Sncg (NM_011430) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Sncg (NM_011430) AAV Particle
Tag: Myc-DDK
Symbol: Sncg
Synonyms: C79089; persyn
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200624 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACGTCTTCAAGAAAGGCTTCTCCATTGCCAAGGAAGGTGTTGTGGGTGCTGTAGAAAAGACCAAGC
AGGGAGTAACGGAGGCAGCTGAGAAGACCAAGGAGGGGTTATGTATGTGGGCACCAAAACCAAGGAGAA
CGTGGTACAAAGTGTACCTCAGTGGCTGAGAAGACCAAGGAGCAGGCCAATGCCGTGAGTGAAGCTGTG
GTCAGCAGCGTCAACACAGTGGCCAACAAGACCGTGGAGGAGGCTGAAAACATCGTGGTACCACCGGGG
TGGTGCCGAAGGAGGACTTGAACCCCTGCACAGGACCAAGAAGCCAAGAGCAAGAGGAGAATGAAGA
GGCCAAGAGTGGAGAAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200624 protein sequence
Red=Cloning site Green=Tags(s)

MDVFKKGFSLAKEGVVGA VEKTKQGVTEAAEKTKEGVMYVGT KTKENVVQSVTSVAEKTKEQANAVSEAV
VSSVNTVANKTVEEAENIVVTTGVVRKEDLEPPAQDQEAKEQEENEAKSGED

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_011430



[View online »](#)

ORF Size:	372 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_011430.1</u>
RefSeq Size:	757 bp
RefSeq ORF:	372 bp
Locus ID:	20618
UniProt ID:	<u>Q9Z0F7</u>
Cytogenetics:	14 20.81 cM
MW:	13.2 kDa