

## Product datasheet for MR200625L4

### BC030500 (NM\_173411) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	BC030500 (NM_173411) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	BC030500
Synonyms:	MGC40983
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200625).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_173411
ORF Size:	372 bp



[View online »](#)

