

## Product datasheet for **MR200631L4V**

### **Bex3 (NM\_001110234) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bex3 (NM_001110234) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bex3
Synonyms:	AL033356; DXWsu67e; Gcap27; Nade; Ngfrap1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001110234
ORF Size:	375 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200631).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001110234.1</a>
RefSeq Size:	893 bp
RefSeq ORF:	345 bp
Locus ID:	12070
UniProt ID:	<a href="#">Q9WTZ9</a>
Cytogenetics:	X 57.64 cM
Gene Summary:	May be a signaling adapter molecule involved in p75NTR-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death.[UniProtKB/Swiss-Prot Function]



[View online »](#)