

Product datasheet for MR200635

Scamp2 (BC003452) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Scamp2 (BC003452) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Scamp2
Synonyms: A1875466; Sc2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR200635 representing BC003452
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGGCCTTCGACACTAACCCCTTCGCGGACCCAGTGGACGTA AACCCCTTCAGGATCCCTCTGTGA
 CCCAGCTGACCAATGCTCCTCAGAGTGGCTGGCTGAGTTCAATCCCTTCTCAGAGACAAATGCAGCGAC
 AACAGTTCCTGCCACACAAGCTCCTGGGCCCTCCAGCCAGCAGTTCTCCAGCCCTCAGTGAACCAGCA
 CAGCCAACGCCCCAGGCTGTTGCAGCTGCGGCCAGGCAGGCTTGCTTCGACAGCATGTGAACTGTACG
 TAGTCGCTGCCGAGCTGGAACGAAGGAGCGCGGTTGCAGCCCTCGGGGGAATTTGCATGTGCGAGA
 TAACGACTGGCCGCACTCCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200635 representing BC003452
 Red=Cloning site Green=Tags(s)

MSAFDTNPFADPVDVNPFDPSVTQLTNAPQSGLAEFNPFSETNAATVPATQAPGPSQPAVLQPSVEPA
 QPTPQAVAAAAQAGLLRQHVELYVVAELERKERGLQPLGGNLHVRDNDWPPLP

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC003452

ORF Size: 372 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC003452</u>
RefSeq Size:	1279 bp
RefSeq ORF:	374 bp
Locus ID:	24044
Cytogenetics:	9 B
MW:	46.9 kDa
Gene Summary:	Functions in post-Golgi recycling pathways. Acts as a recycling carrier to the cell surface (By similarity).[UniProtKB/Swiss-Prot Function]