

Product datasheet for MR200640L3

Qrfp (NM_183424) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Qrfp
Synonyms:	P51; P518; QR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200640).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.

ACCN:	NM_183424
ORF Size:	372 bp



OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_183424.1 , NP_906269.1
RefSeq Size:	2863 bp
RefSeq ORF:	375 bp
Locus ID:	227717
UniProt ID:	Q8CE23
Cytogenetics:	2 B
Gene Summary:	This gene encodes a preproprotein that is proteolytically processed to generate multiple protein products. The encoded products are members of the RFamide family of neuropeptides, characterized by their common protein C-terminus consisting of an arginine (R) and an amidated phenylalanine (F). These products include the neuropeptides QRFP-26 (26RFa) and the N-terminally extended form, QRFP-43 (43RFa). Both of these neuropeptides bind to the pyroglutamylated RFamide peptide receptor (QRFPR) and may regulate blood pressure, reproduction and food intake. [provided by RefSeq, Sep 2015]

Product images:



Circular map for MR200640L3