

Product datasheet for MR200642

Rpl31 (NM_053257) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl31 (NM_053257) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rpl31
Synonyms: Dts
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200642 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCCCGCAAAGAAGGGTGGCGAGAAGAAGGGCCGTTCTGCCATCAACGAGGTGGTGACCCGAG
AATACACCATCAACATTCACAAGCGCATCCATGGAGTGGGCTTCAAGAAGCGTGCTCCTCGGGCACTCAA
AGAAATTCGGAAGTTTGCCATGAAGGAAATGGGGACACCAGATGTGCGCATTGACACCAGGCTCAATAAA
GCCGTCTGGCCAAGGAATAAGGAACGTTCCATATCGCATCCGAGTACGCTTGCCAGAAAACGTAATG
AGGATGAGGATCCCAAACAAGCTCTACACGCTGGTAACCTACGTGCCTGTACCACATTCAAAAATCT
ACAGACGGTCAACGTGGATGAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200642 protein sequence
Red=Cloning site Green=Tags(s)

MAPAKKGGEEKKGRSAINVVREYTIINIKRIHGVGFKKRAPRALKEIRKFAMKEMGTPDVRIDTRLNK
AVWAKGIRNVPYRIRVRLSRKRNEDESPNKLYTLVTVVPTTFKNLQTVNVNEN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_053257

ORF Size: 378 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053257.3](#), [NP_444487.1](#)

RefSeq Size: 1047 bp

RefSeq ORF: 378 bp


Locus ID: 114641

UniProt ID: [P62900](#)

Cytogenetics: 1 B

MW: 14.5 kDa

Product images:



Circular map for MR200642