

## Product datasheet for **MR200648**

### Rps25 (NM\_024266) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rps25 (NM\_024266) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rps25  
**Synonyms:** 2810009D21Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR200648 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCTCCCAAAGACGACAAGAAGAAGAAAGATGCCGAAAGTCGGCCAAAAGGATAAAGACCCAGTAA  
ATAAATCTGGTGGCAAGGCCAAGAAGAAGAGTGGTCCAAAGGCAAAGTTCGGGACAAGTTGAACAATCT  
TGTCTGTTCGACAAAGCGACATACGACAAGCTCTGTAAGGAGGTTCCGAACTATAAGCTTATTACTCCA  
GCCGTGGTCTCTGAGAGACTGAAGATTCGCGGTTCTTGGCCAGGGCAGCCCTTCAGGAGCTCCTTAGTA  
AAGGACTTATCAAGCTGGTTTCAAAGCACAGAGCCCAAGTAATTTACACCAGAAACAAAAGGGTGGGA  
CGCTCCAGCTGCTGGCGAAGATGCA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200648 protein sequence  
Red=Cloning site Green=Tags(s)  
MPPKDDKKKDKAGKSAKKDKDPVKNSSGGKAKKKKWSKGVDRKLNLLVLFDKATYDKLCKEVPNYKLITP  
AVVSERLKIRGSLARAALQELLSKGLIKLVSKHRAQVIYTRNTKGGDAPAAGEDA

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


ACCN: NM\_024266

ORF Size: 378 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_024266.3](#), [NP\\_077228.1](#)

RefSeq Size: 701 bp

RefSeq ORF: 378 bp


Locus ID: 75617

UniProt ID: [P62852](#)

Cytogenetics: 9 A5.2

MW: 13.7 kDa

**Product images:**



Circular map for MR200648