

## Product datasheet for **MR200667**

### Cxcl9 (NM\_008599) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cxcl9 (NM\_008599) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cxcl9  
**Synonyms:** BB139920; CMK; crg-10; Mig; MuMIG; Scyb9  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR200667 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAGTCCGCTGTTCTTTCTTTGGGCATCATCTTCCTGGAGCAGTGTGGAGTTCGAGGAACCCCTAG  
TGATAAGGAATGCACGATGCTCCTGCATCAGCACCAGCCGAGGCACGGTCCACTACAAATCCCTCAAAGA  
CCTCAAACAGTTTGCCCCAAGCCCAATTGCAACAAAACCTGAAATCATTGCTACACTGAAGAACGGAGAT  
CAAACCTGCCTAGATCCGACTCGGCAATGTGAAGAAGCTGATGAAAGAATGGGAAAAGAAGATCAACC  
AAAAGAAAAGCAAAGAGGGGGAAAAACATCAAAGAACATGAAAACAGAAAACCCAAAACACCCCA  
AAGTCGTCGTCGTTCAAGGAAGACTACA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200667 protein sequence  
Red=Cloning site Green=Tags(s)  
MKSAVL FLLGIIFLEQCGV RGLVIRNARCSCISTSRGTVHYKSLKDLKQFAPSPNCNKTEIIATLKNGD  
QTCLDPDSANVKKLMKEWEKKINQKKKQKRGKKHQKMNMRKPKTPQSRRRSRKTT

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_008599

**ORF Size:** 381 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_008599.1](#), [NM\\_008599.2](#), [NM\\_008599.3](#), [NM\\_008599.4](#), [NP\\_032625.2](#)


**RefSeq Size:** 2905 bp

**RefSeq ORF:** 381 bp

**Locus ID:** 17329

UniProt ID: [P18340](#)  
 Cytogenetics: 5 46.51 cM  
 MW: 14.5 kDa  
 Gene Summary: May be a cytokine that affects the growth, movement, or activation state of cells that participate in immune and inflammatory response.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR200667