

Product datasheet for **MR200667L3V**

Cxcl9 (NM_008599) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cxcl9 (NM_008599) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cxcl9
Synonyms:	BB139920; CMK; crg-10; Mig; MuMIG; Scyb9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008599
ORF Size:	381 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200667).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008599.4 , NP_032625.2
RefSeq Size:	2905 bp
RefSeq ORF:	381 bp
Locus ID:	17329
UniProt ID:	P18340
Cytogenetics:	5 46.51 cM
Gene Summary:	May be a cytokine that affects the growth, movement, or activation state of cells that participate in immune and inflammatory response.[UniProtKB/Swiss-Prot Function]



[View online »](#)