

## Product datasheet for **MR200673A1V**

### Mouse Wfdc5 (NM\_145369) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Wfdc5 (NM\_145369) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Wfdc5  
**Synonyms:** BC019734; Prg5; Wap1  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200673 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAGAGTCCAGAGCTCCCTCCTCTGGTAGTCCTCCTGGCTTTGGAGACTCAGTTGCCTGTTGCCTTGT  
 GTAGAAAGAAGGGAGACAAATTGGGGGTTGCCACCAGACGATGGACTCTGCAACCAGAATACTCCTCC  
 TGATCAGTGCTTGAATGACAAACAGTGTCCCTCTTCTTGAAGTGCTGTGCGCGAGCCTGCTTCTCCAG  
 TGTATGCCTCGTGTCTTTGTAAGTCGGGCAAATGCCCGTGGATCAACTTCGCTGCCTCAGTCCCACAA  
 AGCACATGTGTAACAAAGACTCGGACTGCTCAGGCAAGAAGCGGTGCTGTGCCAGCGCCTGTGGCCGGGA  
 CTGCCGAGACCCAGCAAAGTACCATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200673 protein sequence  
 Red=Cloning site Green=Tags(s)

MRVQSSLLLLVLLALETQLPVALCRKKGDKLGGCPPDDGLCNQNTPPDQCLNDKQCPSSWKCCRACFLQ  
 CMPRVFVKSGKCPVDQLRCLSPKHMCKNDSDCSGKKRCCASACGRDCRDPSTI

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_145369



[View online »](#)

<b>ORF Size:</b>	381 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_145369.2</a></u>
<b>RefSeq Size:</b>	2324 bp
<b>RefSeq ORF:</b>	381 bp
<b>Locus ID:</b>	209232
<b>Cytogenetics:</b>	2 H3
<b>MW:</b>	13.8 kDa