

## Product datasheet for **MR200677A1V**

### Mouse Nutf2 (NM\_026532) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Nutf2 (NM\_026532) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Nutf2  
**Synonyms:** 2700067I02Rik; AI115459; AI413596; AW546000; Ntf2; Ntf3; PP15  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200677 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGAGACAAGCCAATTTGGGAGCAGATTGGATCCAGCTTTATTCAACATTACTACCAGTTATTTGATA  
ACGACAGAACCCAAGTAGGCGCAATTTATTTGATGCATCATGCCTTACATGGGAAGGACAGCAGTTCCA  
GGGAAAGCTGCCATTGTGGAGAAGCTATCTAGCCTCCGTTCCAGAAAATTAGCATAGCATCACGGCG  
CAGGACCACAGCCTACACCAGATAGCTGCATCATCAGCATGGTTGTGGCCAGCTCAAGGCTGATGAAG  
ATCCCATCATGGGTTTCCACCAGATGTTTCTATTAAGAACATCAACGATGCTTGGGTTTGCACCAATGA  
CATGTTCCAGGCTTGCCCTGCACAACCTTCGGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200677 protein sequence  
 Red=Cloning site Green=Tags(s)

MGDKPIWEQIGSSFIQHYYQLFDNDRTQLGAIYIDASCLTWEGQQFQGKAAIVEKLSLPPFKIQHSITA  
QDHQPTDSCIISMVVGQLKADEDPIMGFHQMFLKKNINDAWVCTNDMFRLALHNF

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_026532



[View online »](#)

<b>ORF Size:</b>	384 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_026532.3</a></u>
<b>RefSeq Size:</b>	1966 bp
<b>RefSeq ORF:</b>	384 bp
<b>Locus ID:</b>	68051
<b>UniProt ID:</b>	<u><a href="#">P61971</a></u>
<b>Cytogenetics:</b>	8 53.05 cM
<b>MW:</b>	14.5 kDa