

Product datasheet for **MR200679**

3830406C13Rik (NM_178141) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 3830406C13Rik (NM_178141) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 3830406C13Rik
Synonyms: 3110059I09Rik; AV167947
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200679 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTTCCTTGTACTCAAGAAATTCACCTTTCTAAAAGACATGAAGAAATATTATCACAAAGATTAA
TGTTACTTCAGAAAATGAAAACAAATTTGGAGATGAAAACACAGAAAGGCATCTCTCTCCAAGCCAC
TGAGACTGCTTCTAGAAGAACTCCGACTCTTGAAGGACATAGATGCTGCAGAGGAGGCTTTTCAGACC
AAGCTTATCCACACCCACAGCCTTCGATGCTTTCTCTTGAGACTCGTTACTGGGCATCAGTAGAAGAAC
ATATCCCAAGTGGAACTGTTTCTTTTAGGAAGAGCACCATATCCAATTGGTGCTGAAAATCAAATGA
GGTACCATTCGTTCAAACCGAGGCAAAGTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200679 protein sequence
Red=Cloning site Green=Tags(s)
MASLFTQEIHLSCRHEEILSQLMLLQKMKNFQDENTERASLLQATETASRRNLRLKIDIAEEAFQT
KLIPHPQPSMLSLETRYWASVEEHIPKWELFLGRAPYPIGAENQNEVPFVQTEAKL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_178141

ORF Size: 384 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178141.3](#)

RefSeq Size: 2422 bp

RefSeq ORF: 384 bp


Locus ID: 218734

UniProt ID: Q8BGD0

Cytogenetics: 14 A1

MW: 14.7 kDa

Product images:



Circular map for MR200679