

Product datasheet for MR200684L3

Mrap (NM_029844) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrap (NM_029844) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Mrap
Synonyms:	1110025G12Rik; C21ORF61; Falp; ORF61
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200684).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.

ACCN:	NM_029844
ORF Size:	381 bp



OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_029844.1
RefSeq Size:	992 bp
RefSeq ORF:	384 bp
Locus ID:	77037
UniProt ID:	Q9D159
Cytogenetics:	16 C3.3
Gene Summary:	Modulator of melanocortin receptors (MC1R, MC2R, MC3R, MC4R and MC5R). Acts by increasing ligand-sensitivity of melanocortin receptors and enhancing generation of cAMP by the receptors. Required both for MC2R trafficking to the cell surface of adrenal cells and for signaling in response to corticotropin (ACTH). May be involved in the intracellular trafficking pathways in adipocyte cells (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR200684L3