

## Product datasheet for **MR200684L3V**

### Mrap (NM\_029844) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mrap (NM_029844) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mrap
Synonyms:	1110025G12Rik; C21ORF61; Falp; ORF61
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029844
ORF Size:	384 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200684).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_029844.1</a>
RefSeq Size:	992 bp
RefSeq ORF:	384 bp
Locus ID:	77037
UniProt ID:	<a href="#">Q9D159</a>
Cytogenetics:	16 C3.3


[View online »](#)

**Gene Summary:**

Modulator of melanocortin receptors (MC1R, MC2R, MC3R, MC4R and MC5R). Acts by increasing ligand-sensitivity of melanocortin receptors and enhancing generation of cAMP by the receptors. Required both for MC2R trafficking to the cell surface of adrenal cells and for signaling in response to corticotropin (ACTH). May be involved in the intracellular trafficking pathways in adipocyte cells (By similarity).[UniProtKB/Swiss-Prot Function]