

Product datasheet for **MR200705A1V**

Mouse Etfbkmt (BC046911) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Etfbkmt (BC046911) AAV Particle
Tag: Myc-DDK
Symbol: Etfbkmt
Synonyms: 4833442J19Rik; AI256744; ETFB-KMT
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200705 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGTGGAGCATCAAAGATCTTGGCCAACGACATAGACCCCATTCAGGAATGGCGATTACCTAAATT
 GCAAATTAATGGACTGAATCCTTTTCCTGTTTTAACTAAAAACATTTTGAATACCAACAAGGTAAGTT
 TGATCTCATTGTTCTTGGAGATATGTTTTATGATGAAGACCTTGACAGACAGTCTTCATCTGTGGCTGCAG
 AACTACTTTTGGACCCATGGAACCCGAGTACTGATTGGTGACCCTGGGAGGCCCCAGTTCAGTGGACATA
 GCATTCGGCATCAACTGTACCAACTAGTAGAATACACACTTCCAGAGCCAACCTCAGCAAGAGAACAATGG
 ACTGACCACAAGTGCTGTATGGGATTTTCATCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200705 protein sequence
 Red=Cloning site Green=Tags(s)

MSGASKILANDIDPIAGMAITLNCKLNGLNPFVLTKNILNTQQGKFDLIVLGDIFYDEDLADSLHLWLQ
 NYFWTHGTRVLIGDPGRPQFSGHSIRHQLYQLVEYTLPEPTQQENGLTTSVWDFHP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC046911



[View online »](#)

ORF Size:	384 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC046911</u> , <u>AAH46911</u>
RefSeq Size:	1104 bp
RefSeq ORF:	386 bp
Locus ID:	320204
Cytogenetics:	6 G3
MW:	14.4 kDa