

## Product datasheet for **MR200713A1V**

### Mouse Hypk (NM\_026318) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Hypk (NM\_026318) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Hypk  
**Synonyms:** 2310003F16Rik; HSPC136  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200713 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGCGCGCGGGTGGAGATCGAAATGGCGACCGAGGGAGACGTGGAGTTGGAGTTAGAGACCGAGACCA  
GCGGCCCGAGCGGCCTCCCGAGAAGCCAAGGAAACATGACAGTGGTGTCTGCTGACCTAGAGCGGGTTAC  
TGACTATGCGGAGGAGAAGGAGATCCAGAGTTCGAATCTTGAACGGCTATGTCCGTATTGGAGACAGA  
CGGTCCAGGGAACAAAAGGCCAAACAAGAGCGGGAGAAGGAAGTGGCGAAGGTCACTATCAAGAAGGAAG  
ATCTGGAGTTGATAATGACGAAATGGAGATCTCTCGAGCAGCAGAAAGGAGCTTGAGGGAACACAT  
GGCAATGTGGTGGAGGCTCTTATTGCCCTAACCAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200713 protein sequence  
Red=Cloning site Green=Tags(s)

MRRRGEIEMATEGDVELELETETSGPERPPEKPRKHDSGAADLERVTDYAEKEIQSSNLETAMSVIGDR  
RSREQAKQEREKELAKVTIKKEDLELIMTEMEISRAAAERSLREHMGNVVEALIALTN

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_026318



[View online »](#)

|                      |   |
|----------------------|---|
| <b>ORF Size:</b>     | 390 bp  |
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68  |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <u><a href="#">NM_026318.1</a></u>  |
| <b>RefSeq Size:</b>  | 620 bp  |
| <b>RefSeq ORF:</b>   | 390 bp  |
| <b>Locus ID:</b>     | 67693   |
| <b>UniProt ID:</b>   | <u><a href="#">Q9CR41</a></u>   |
| <b>Cytogenetics:</b> | 2 E5  |
| <b>MW:</b>           | 14.7 kDa  |