

## Product datasheet for **MR200717L4V**

### 6230409E13Rik (BC059864) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	6230409E13Rik (BC059864) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	6230409E13Rik
Synonyms:	6230409E13Rik; AW742931
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC059864
ORF Size:	387 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200717).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC059864</a> , <a href="#">AAH59864</a>
RefSeq Size:	4035 bp
RefSeq ORF:	389 bp
Locus ID:	76132
Cytogenetics:	4 A3
Gene Summary:	May play a role in axonal development.[UniProtKB/Swiss-Prot Function]


[View online »](#)