

Product datasheet for **MR200724L4V**

Aim2 (NM_001013779) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aim2 (NM_001013779) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aim2
Synonyms:	Gm1313; lfi210
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001013779
ORF Size:	390 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200724).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013779.1
RefSeq Size:	2541 bp
RefSeq ORF:	1065 bp
Locus ID:	383619
UniProt ID:	Q91VJ1
Cytogenetics:	1 H3



[View online »](#)

Gene Summary:

Involved in innate immune response by recognizing cytosolic double-stranded DNA and inducing caspase-1-activating inflammasome formation in macrophages. Upon binding to DNA is thought to undergo oligomerization and to associate with PYCARD initiating the recruitment of caspase-1 precursor and processing of interleukin-1 beta and interleukin-18. Detects cytosolic dsDNA of viral and bacterial origin in a non-sequence-specific manner. Can also trigger PYCARD-dependent, caspase-1-independent cell death that involves caspase-8. [UniProtKB/Swiss-Prot Function]