

## Product datasheet for **MR200733A1V**

### Mouse Cdkn2b (NM\_007670) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Cdkn2b (NM\_007670) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Cdkn2b  
**Synonyms:** AV083695; INK4b; MTS2; p15; p15(INK4b); p15INK4b  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200733 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTTGGGCGGCAGCAGTGACGCGGGCCTGGCCACCGCCGCGCGGGGCAAGTGGAGACGGTGCGGC  
AGCTCCTGGAAGCCGGCGCAGATCCCAACGCCCTGAACCGCTTCGGGAGGCGCCCAATCCAGGTCATGAT  
GATGGGCGAGCGCCAGGTGGCAGAGCTGCTGCTGCTCCACGGAGCAGAACCCTGCGCCGACCTGCC  
ACCTTACCAGACCTGTGCACGACGAGCTCGGGAAGGCTTCTGGACACGCTTGTGCTGCACCGGG  
CAGGGGCGGGCTGGATGTGTGACGCTCGGGCCGCTCGCGGTAGACTTAGCTGAAGAGCAGGGCCA  
CCGTGACATTGCGAGGTATCTGCACGCTGCCACTGGAGAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200733 protein sequence  
Red=Cloning site Green=Tags(s)

MLGSSDAGLATAAARGQVETVRQLLEAGADPNALNRFGRRIQVMMGSAQVAELLLLHGAEPNCADPA  
TLTRPVHDAAREGFLDTLVVLRHAGARLDVCDAWGRLPVDLAEEQGHRDIARYLHAATGD

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_007670



[View online »](#)

<b>ORF Size:</b>	393 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_007670.3</a></u>
<b>RefSeq Size:</b>	1393 bp
<b>RefSeq ORF:</b>	393 bp
<b>Locus ID:</b>	12579
<b>UniProt ID:</b>	<u><a href="#">P55271</a></u>
<b>Cytogenetics:</b>	4 42.15 cM
<b>MW:</b>	13.8 kDa