

Product datasheet for **MR200737A1V**

Mouse 1700012A03Rik (NM_029587) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse 1700012A03Rik (NM_029587) AAV Particle
Tag: Myc-DDK
Symbol: 1700012A03Rik
Synonyms: AI429084
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200737 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATACTGACATGAGACGAGCATACACTGTTATTCATTCACGTGTCCATGCCTCCTGGCAAGGGTGCC
TTGACCTTCAGGAGCTTCCACTGAGTTCAAGATCTGGAGTCACTCTAGAGGAGCGGTGAGTAAGAGAAA
GCAAAGAGTCTGCAGAAACAAGCCTGCTGCAGAAGAGGGCGCAGAGAGAAGCCTGCTGCAGAAGTGGCTG
CAGAGAGATGCCTGCTGCAGAGAGAAGCCTCCTGGCCAGAGAAGCCTGCTGGAGAAGTGGCTGCAGAGAG
ATGCCTGCTGCAGAAATGACTGCAGAGAGGAGAAGGTTGCTGTAAAAGCCACTGTAGAGACAAGTGTAGC
ACCTTATCAGGGCCAAGAAGCCACCTTCTCACCTGACCCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200737 protein sequence
Red=Cloning site Green=Tags(s)

MDTDMRRAYTVIHSRVHASWQGCVDLQELSTEFKIWSHSRGAVSKRKQRVCRNKPAEEGAERSLLQKWL
QRDACCREKPPGQRSLLLEKWLQRDACCRRNDCEEKVAVKATVETSVAPYQGQEATFSPDP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_029587



[View online »](#)

ORF Size:	393 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_029587.1</u>
RefSeq Size:	918 bp
RefSeq ORF:	393 bp
Locus ID:	76382
Cytogenetics:	6 A3.3
MW:	14.9 kDa