

Product datasheet for **MR200742A1V**

Mouse Dnajc19 (BC024335) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Mouse Dnajc19 (BC024335) AAV Particle
 Tag: Myc-DDK
 Symbol: Dnajc19
 Synonyms: 1810055D05Rik; AA959924; Tim14
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >MR200742 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCCAGCACAGTGGTAGCAGTCGGGTTGACCATTGCTGCTGCAGGATTTGCAGGCCGTTATGTTTTAC
 AAGCCATGAAGCATGTGGAGCCTCAAGTAAACAAGTTTTTCAGAGCCTACCAAATCTGCATTCGGTGG
 TGGGTACTACAGAGGTGGATTTGAACCCAAATGACAAACGGGAAGCAGCATTAAATATTAGGTGAAGC
 CCTACTGCCAATAAAGGAAGATCAGGGATGCTCATCGCCGATTATGCTATTAAATACCCAGACAAGG
 GTAAGCAACTTTTATTACTATACTGGACTGTTAATGGAATAGCTAATTCTTGGTGTATGAATGTCTGTG
 TGTGAGCAGTTTCCAGTTGTGGATCACATTGTAGTGTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200742 protein sequence
 Red=Cloning site Green=Tags(s)

MASTVVAVGLTIAAAGFAGRYVLQAMKHVEPQVKQVFQSLPKSAFGGGYRGGFEPKMTKREAALILGVS
 PTANKGKIRDAHRRIMLLNHPDKGKQLLLLYWTVNGIANSWCYECLCVSSFPPVDHIVVS

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
 Serotype: AAV-2
 ACCN: BC024335



[View online »](#)

ORF Size:	390 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC024335 , AAH24335
RefSeq Size:	1716 bp
RefSeq ORF:	392 bp
Locus ID:	67713
Cytogenetics:	3
MW:	14.2 kDa