

Product datasheet for **MR200747L3V**

Ndufa6 (NM_025987) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ndufa6 (NM_025987) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ndufa6
Synonyms:	14kDa; 2700038D15Rik; B230217P19Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025987
ORF Size:	393 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200747).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025987.3 , NP_080263.1
RefSeq Size:	560 bp
RefSeq ORF:	396 bp
Locus ID:	67130
UniProt ID:	Q9CQZ5
Cytogenetics:	15 E1



[View online »](#)

Gene Summary:

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed to be not involved in catalysis. Required for proper complex I assembly. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone. [UniProtKB/Swiss-Prot Function]