

## Product datasheet for **MR200750L4V**

### Slpi (NM\_011414) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slpi (NM_011414) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slpi
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011414
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200750).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011414.1</a>
RefSeq Size:	894 bp
RefSeq ORF:	396 bp
Locus ID:	20568
UniProt ID:	<a href="#">P97430</a>
Cytogenetics:	2 H3



[View online »](#)

**Gene Summary:**

Acid-stable proteinase inhibitor with strong affinities for trypsin, chymotrypsin, elastase, and cathepsin G (PubMed:9126337). Modulates the innate immune response after bacterial infection (PubMed:12615907). Contributes to regulate the inflammatory and immune responses to the intracellular parasite *L.major* (PubMed:25030421). Down-regulates responses to bacterial lipopolysaccharide (LPS) (PubMed:9039268, PubMed:12615907, PubMed:25030421). Plays a role in regulating the activation of NF-kappa-B and inflammatory responses (PubMed:11017147, PubMed:12615907). Has antimicrobial activity against mycobacteria, but not against salmonella (PubMed:18322212). Contributes to normal resistance against infection by *M.tuberculosis* (PubMed:18322212). Required for normal resistance to *L.major* (PubMed:25030421). Required for normal wound healing, probably by preventing tissue damage by limiting protease activity (PubMed:11017147, PubMed:25030421). Together with ELANE, required for normal differentiation and proliferation of bone marrow myeloid cells (By similarity).[UniProtKB/Swiss-Prot Function]