

Product datasheet for **MR200754L3V**

Mmgt1 (NM_146234) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mmgt1 (NM_146234) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mmgt1
Synonyms:	9630048L06Rik; BC032271; Emc5; mir-5116; Mir5116; Tmem32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_146234
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200754).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_146234.2
RefSeq Size:	4084 bp
RefSeq ORF:	396 bp
Locus ID:	236792
UniProt ID:	Q8K273
Cytogenetics:	X A6
Gene Summary:	Mediates Mg(2+) transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)