

Product datasheet for **MR200762A1V**

Mouse Fabp4 (NM_024406) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Fabp4 (NM_024406) AAV Particle
Tag: Myc-DDK
Symbol: Fabp4
Synonyms: 422/aP2; ALBP/Ap2; Ap2; Lbpl
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200762 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGTGATGCCTTTGTGGGAACCTGGAAGCTTGTCTCCAGTGAAAACCTCGATGATTACATGAAAGAAG
TGGGAGTGGGCTTTGCCACAAGGAAAGTGGCAGGCATGGCCAAGCCCAACATGATCATCAGCGTAAATGG
GGATTTGGTCACCATCCGGTCAGAGAGTACTTTAAAAACACCGAGATTCCTTCAAACCTGGGCGTGAA
TTCGATGAAATCACCGCAGACGACAGGAAGGTGAAGAGCATCATAACCCTAGATGGCGGGGCCCTGGTGC
AGGTGCAGAAGTGGGATGGAAAGTCGACCACAATAAAGAGAAAACGAGATGGTGACAAGCTGGTGGTGA
ATGTGTTATGAAAGCGTGACTTCCACAAGAGTTTATGAAAGGGCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200762 protein sequence
Red=Cloning site Green=Tags(s)

MCDAFVGTWKLVSSENFDDYMKEVGVGFATRKYVAGMAKPNMIIISVNGDLVTIRSESTFKNTEISFKLGVE
FDEITADDRKVKSIITLDGGALVQVQKWDGKSTTIKRKRDKLVEECVMKGVSTRVYERA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_024406



[View online »](#)

ORF Size:	399 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024406.1</u>
RefSeq Size:	646 bp
RefSeq ORF:	399 bp
Locus ID:	11770
UniProt ID:	<u>P04117</u>
Cytogenetics:	3 2.56 cM
MW:	14.6 kDa