

## Product datasheet for **MR200775L4V**

### Lsm1 (NM\_026032) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lsm1 (NM_026032) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lsm1
Synonyms:	2810025O06Rik; CASM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026032
ORF Size:	402 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200775).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026032.1</a> , <a href="#">NP_080308.1</a>
RefSeq Size:	2360 bp
RefSeq ORF:	402 bp
Locus ID:	67207
UniProt ID:	<a href="#">Q8VC85</a>
Cytogenetics:	8 A2



[View online »](#)

**Gene Summary:**

Plays a role in the degradation of histone mRNAs, the only eukaryotic mRNAs that are not polyadenylated (By similarity). Probably also part of an LSM subunits-containing complex involved in the general process of mRNA degradation (By similarity).[UniProtKB/Swiss-Prot Function]