

Product datasheet for **MR200781L3V**

Rps24 (BC091748) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rps24 (BC091748) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rps24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC091748
ORF Size:	399 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200781).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC091748 , AAH91748
RefSeq Size:	522 bp
RefSeq ORF:	401 bp
Locus ID:	20088
Cytogenetics:	14 A3
Gene Summary:	Required for processing of pre-rRNA and maturation of 40S ribosomal subunits. [UniProtKB/Swiss-Prot Function]



[View online »](#)