

Product datasheet for **MR200794**

Cck (BC028487) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cck (BC028487) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cck
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200794 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAGCGGCGTATGTCTGTGCGTGGTATGGCAGTCCTAGCTGCTGGCGCCCTGGCGCAGCCGGTAG
TCCCTGCAGAAGCTACGGACCCCGTGGAGCAGCGGGCGCAAGAGGCGCCCGAAGGCAGCTGCGGGCTGT
GCTCCGGACGGACGGCGAGCCCCGAGCGCGCTGGCGCACTGCTAGCGGATACATCCAGCAGGTCCGC
AAAGTGGCATGGATGGTACCTCTGGTTGGGTGTTAACCTGGACCAGTAGAGCTGGATTAACACAGGA
GGTGGGCATCCTTTCTGTGGAGCTCACGAACCAATTTTCTGCCCGCATTTGAACAGCCCATGGCTTG
TCGACCTGTCTGCATTTGGCTTGACTGCTCCTTCTGGCCGCATGTCCGTTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200794 protein sequence
Red=Cloning site Green=Tags(s)

MKSGVCLCVVMAVLAAGALAQPVVPAEATDPVEQRAQEAPRRQLRAVLRDGEPRARLGALLARYIQQVR
KVAWMVTSGWVLTWTSRAGLKHRRWASFLWSSRTQFFLPAFEQPMACRPVCIWLDCSFWPHVRS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: BC028487

ORF Size: 402 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC028487](#), [AAH28487](#)

RefSeq Size: 816 bp

RefSeq ORF: 404 bp


Locus ID: 12424

Cytogenetics: 9 F4

MW: 15.2 kDa

Gene Summary: This gene encodes a member of the gastrin/cholecystokinin family of proteins. The encoded preproprotein is proteolytically processed to generate multiple protein products, including the peptide hormones cholecystokinin-8, -12, -33, and others. The encoded peptides have been shown to regulate gastric acid secretion and food intake. A sulfated form of cholecystokinin-8 may modulate neuronal activity in the brain. Homozygous knockout mice for this gene exhibit impaired insulin secretion, enhanced insulin sensitivity, and resistance to obesity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2015]

Product images:



Circular map for MR200794