

Product datasheet for **MR200794A1V**

Mouse Cck (BC028487) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Cck (BC028487) AAV Particle
Tag: Myc-DDK
Symbol: Cck
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200794 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAGCGCGTATGTCTGTGCGTGGTGATGGCAGTCCTAGCTGCTGGCGCCCTGGCGCAGCCGGTAG
TCCCTGCAGAAGCTACGGACCCCGTGGAGCAGCGGGCGCAAGAGGCGCCCGAAGGCAGCTGCGGGCTGT
GCTCCGGACGGACGGCGAGCCCCGAGCGCGCTGGGCGCACTGCTAGCGGATACATCCAGCAGGTCCGC
AAAGTGGCATGGATGGTACCTCTGGTGGGTGTTAACCTGGACCAGTAGAGCTGGATTAACACACAGGA
GGTGGGATCCTTTCTGTGGAGCTCACGAACCAATTTTCTGCCCGCATTTGAACAGCCCATGGCTTG
TCGACCTGTGCATTTGGCTTGACTGCTCCTTCTGGCCGCATGTCCGTTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200794 protein sequence
Red=Cloning site Green=Tags(s)

MKSGVCLCVMAVLAAGALAQPVVPAEATDPVEQRAQEAPRRQLRAVLRDGEPRARLGALLARYIQQVR
KVAWMVTSGWLTWTSRAGLKHRRWASFLWSSRTQFFLPAFEQPMACRPVCIWLDCSFVPHVRS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC028487
ORF Size: 402 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC028487</u> , <u>AAH28487</u>
RefSeq Size:	816 bp
RefSeq ORF:	404 bp
Locus ID:	12424
Cytogenetics:	9 F4
MW:	15.2 kDa