

## Product datasheet for **MR200798**

### Mturn (NM\_026629) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mturn (NM\_026629) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Mturn  
**Synonyms:** 2410066E13Rik; AI448550; B230212L03Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR200798 representing NM\_026629  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCGGGAGAAGTGGTGTCCAGCACGCCCTTTGAGCTCATCGCCGCCGAGGAGACGGAGCGCAGGAT  
GGACTTCTACGCCGACCCGGGCGTCTCTTCTACGTACTGTGTCCGGACAACGGCTGCGGCGACAGCTTT  
CACGTGTGGAGCGAGAGCGAGGACTGCCTGCCCTTCTGCAGCTGGCCCAAGATTACATCTCTTCTGCG  
GCAAGAAGACCCTCCATGAAGTCTTGAAAAAGTCTTCAAATCCTTCAGACCTTTGCTGGGGCTTCCAGA  
CGCAGATGACGACGCGTTTGAAGAGTACAGTCCGATGTGGAGGAGGAGCCAGAGCGGACCCCG  
CAAATGGGGTCACTCAGCAGTAACTCCTAGGTACCTTGCAGAAGGGGAG


AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9051\\_g04.zip](https://cdn.origene.com/chromatograms/mm9051_g04.zip)

**Restriction Sites:** SgfI-RsrII



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_026629

**ORF Size:** 402 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_026629.3, NP\\_080905.1](#)  
 RefSeq Size: 1710 bp  
 RefSeq ORF: 405 bp  
 Locus ID: 68235  
 Cytogenetics: 6 B3  
 MW: 15.3 kDa  
 Gene Summary: May be involved in early neuronal development.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR200798