

Product datasheet for **MR200801A1V**

Mouse Tmem100 (NM_026433) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Tmem100 (NM_026433) AAV Particle
Tag: Myc-DDK
Symbol: Tmem100
Synonyms: 1810057C19Rik; AV011897
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200801 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACCGAAGAATCCACAAAAGAGAACCTGGGAGCTCCAAAATCTCCACACCTGTGACAATGGAGAAAA
 ACCCCAAGAGGGAAGTTGTGGTCACCACGGGACCCTTGGTCAGCGAGGTTGAGCTGATGGCCGCCACCGG
 GGGTGCCGAACCTCCTGCTACCGCTGCATCATCCCTTTGCCGTGGTGGTCTTCATCACTGGGATTGTG
 GTCACCGCTGTGGTTACAGCTTCAATTCCCATGGTTCCATCATCTCCATCTTCGGCCTGGTCCCTTCTGT
 CCTCCGGACTGTTTTACTAGCCTCCAGTGCCTTGTGCTGGAAGGTGAGACAAAGGAACAAGAAAGTCAA
 GAGACGCGAGAGTCAGACCGCTCTTGTGGTAAATCAGAGGTGCTTGTGGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200801 protein sequence
 Red=Cloning site Green=Tags(s)

MTEESTKENLGAPKSPTPVTMEKNPKREVVVTTGPLVSEVQLMAATGGAELSCYRCIIPFAVVVFITGIV
 VTAVAYSFNHSGSIISIFGLVLLSSGLFLLASSALCWKVRQRNKKVKRRESQ TALVNVNQRCLFA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_026433



[View online »](#)

ORF Size:	405 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026433.1</u>
RefSeq Size:	1764 bp
RefSeq ORF:	405 bp
Locus ID:	67888
UniProt ID:	<u>Q9CQG9</u>
Cytogenetics:	11 C
MW:	14.5 kDa