

## Product datasheet for **MR200802A1V**

### Mouse Rps17 (NM\_009092) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Rps17 (NM\_009092) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Rps17  
**Synonyms:** MGC107578; MGC107579  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200802 representing NM\_009092  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGCCCGGTTTCGCACCAAGACCGTGAAGAAGGCTGCCCGGTCATCATCGAGAAGTACTACACGCGTC  
TGGGTAATGACTTCCACACCAACAAGCGCGTGTGCGAGGAGATCGCCATTATCCCAGCAAGAAGCTTCC  
GAACAAGATAGCCGGCTATGTCACGCATCTGATGAAGCGGATTCAGAGAGGGCCTGTGAGAGGCATCTCT  
ATCAAGTTGCAGGAAGAAGAGAGAGAGAGAGAGATAACTACGTTCTGAGGTCTCAGCCCTAGATCAGG  
AGATCATTGAGGTGGATCCCACACCAAGGAGATGCTGAAGCTCCTGGACTTTGGCAGTCTCTAACCT  
TCAGGTCACTCAGCCTACAGTTGGCATGAATTTCAAAACACCACGTGGAGCTGTT

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200802 representing NM\_009092  
Red=Cloning site Green=Tags(s)

MGRVTRTKVKKAARVIEKYYTRLGNDFHTNKRVCEEIAIIPSKLRNKIAGYVTHLMKRIQRGPVVRGIS  
IKLQEEERERRDNYVPEVSALDQEIIIEVDPDTKEMLKLLDFGSLSNLQVTQPTVGMNFKTPRGAV

**SGP**TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_009092



[View online »](#)

<b>ORF Size:</b>	405 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_009092.3</a>, <a href="#">NP_033118.1</a></u>
<b>RefSeq Size:</b>	489 bp
<b>RefSeq ORF:</b>	408 bp
<b>Locus ID:</b>	20068
<b>UniProt ID:</b>	<u><a href="#">P63276</a></u>
<b>Cytogenetics:</b>	7 D3
<b>MW:</b>	16 kDa