

Product datasheet for MR200825

Rpl27 (NM_011289) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl27 (NM_011289) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rpl27
Synonyms: MGC103134
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200825 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGGGCAAGTTCATGAAACCCGGGAAAGTGGTCTGGCTGGCCGGACGCTACTCCGGACGCAAAGCCG
 TCATCGTGAAGAACATTGACGATGGCACCTCGGACCGCCCTTACAGCCATGCCCTGGTGGCTGGAATTGA
 CCGCTATCCCCGGAAAGTGACAGCTGCCATGGGCAAGAAGAAGATCGCCAAGCGATCCAAGATCAAGTCC
 TTTGTGAAAGTTTATAACTACAACCACCTCATGCCACAAGTACTCTGTGGATATCCCCTTGGACAAGA
 CTGTTGTCAACAAGGATGTGTTCAAGGACCCAGCTTTGAAGCGCAAGGCCAGGCGGGAGGCCAAGGTCAA
 GTTTGAGGAGCGATACAAGACAGGGAAGAACAATGGTTTTCCAGAAGCTTCGCTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200825 protein sequence
 Red=Cloning site Green=Tags(s)

MGKFMKPGKVVLVLAGRYSGRKAVIVKNIDGTSRDPYSHALVAGIDRYPRKVTAMGKKKIAKRSKIKS
 FVKVYNYNHLMPTRYSDIPLDKTVVNKDVFDRPALKRKARREAVKFEERYKTGKNKWFQKLR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_011289

ORF Size: 411 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011289.3](#)


RefSeq Size: 734 bp

RefSeq ORF: 411 bp

Locus ID: 19942

UniProt ID: [P61358](#)
 Cytogenetics: 11 D
 MW: 15.8 kDa
 Gene Summary: Component of the large ribosomal subunit (By similarity). Required for proper rRNA processing and maturation of 28S and 5.8S rRNAs (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR200825