

Product datasheet for **MR200831A1V**

Mouse Nradd (BC039808) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Nradd (BC039808) AAV Particle
Tag: Myc-DDK
Symbol: Nradd
Synonyms: 2610311B09Rik; BB104149; NRDD; NRH2; PLAIDD
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200831 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTTTATAACGTCAGCAAAGGTGTGGTCTATTCAGATACAGCCCTGCAGGGGCAGGACGGGGACAGGG
AAGGAATGTGGGTAGGAGCTGGGGGAGCCCTAGCCCCAATACCTCCTCCCTATTTCCCCCTGAGCCTCC
AGGGCCCTCGAGCAACATCATTCTGTCTACTGTCTCCTAGCTACAGTGATCCTTGGTCTGTGGCC
TATGTGGCCTTCAAATGCTGGCGCTCACATAAGCAAAGGCAACAGTTGGCTAAAGCTCGGACTGTAGAGC
TAGGGGACCCTGACAGGGACCAGAGGCGTGGTGACAGCAACGTCTTCGTGGACTCTCCTCTAGTCTGGA
GCCCTGTATCCAGCCAGGGACCACACCCGGACCTTGGCTGCCAGCCTATACCCCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200831 protein sequence
Red=Cloning site Green=Tags(s)

MLYNVSKGVVYSDTALQGQDGDREGMWVGAGGALAPNTSSLFPPEPPGASSNIIPVYCALLATVILGLLA
YVAFKCRWSHKQRQQLAKARTVELGDPDRDQRRGDSNVFVDSPPSLEPCIPSQGPLDLGCQPIPC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC039808



[View online »](#)

ORF Size:	408 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC039808</u> , <u>AAH39808</u>
RefSeq Size:	1059 bp
RefSeq ORF:	410 bp
Locus ID:	67169
Cytogenetics:	9 F2
MW:	14.5 kDa