

Product datasheet for MR200832L4

Sft2d1 (BC006777) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sft2d1 (BC006777) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Sft2d1
Synonyms:	5630401J11Rik; AA409698
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200832).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	BC006777
ORF Size:	408 bp



[View online »](#)

