

## Product datasheet for MR200839

### 1700021F07Rik (NM\_028158) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 1700021F07Rik (NM\_028158) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** 1700021F07Rik  
**Synonyms:** RP23-18F9.1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR200839 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGCAGAAACCGCTCAGCACCGGCCGAGAGCGCATGAACCTTGTGGCCAGGATGAGATTTGGA  
 AATACCGTCTGAGGGCCGAAAGTGAAGCCCGCAAACCTGGCCGGCAAGTGGGGATATTTAACACGTC  
 CATGAAAGAGCTGCTTGAGGGTGAAGAGGAGCCACAAACCCGAAGCCCAAGCCTGAGCTGCCTAGCCAC  
 TTCTATGTCCGTCAGTGAGCCCATGGATAAACACATCAAGATTCTCCCATCGCCCCAGTCCCAAGA  
 CCACCCAGGGCTTCATCGGATGGCGATCCGGGAAGCCAGCCCTGTACTGTTGGAGAAGTACTCAGAGGT  
 TTGCAGCTGCAAAGGGGCTTATGCCGAGAGCTATGCTGGCCGAGCAAGGCGTGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200839 protein sequence  
 Red=Cloning site Green=Tags(s)

MAQKPLSTTAAERMNLVAQDEIWKYRLRAESEARQNWPAKWGYLTTSMKELLEGEETPQTPKPKPELPSH  
 FYYRVPSPMDKHILPSPVPKTTQGF IGWRSKGPALYCLEKYSEVCSCKGAYARELCWPEQGVH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_028158

**ORF Size:** 411 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_028158.1](#), [NP\\_082434.1](#)

**RefSeq Size:** 641 bp

**RefSeq ORF:** 411 bp


**Locus ID:** 72221

UniProt ID: Q9DA42

Cytogenetics: 2 H4

MW: 15.6 kDa

**Product images:**



Circular map for MR200839