

Product datasheet for MR200842

Ndufb7 (NM_025843) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ndufb7 (NM_025843) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ndufb7
Synonyms: 1110002H15Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200842 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGGCGCACCTGACCCGGCGCTATCTGTGGGATGCCTCGGTGGAGCCTGACCCGAGAAGATACCCA
 GCTTCCC GCCGACCTCGGCTTCCGGAGCGCAAGGAGCGAGTGATGGTGGCCACACAACAAGAGATGAT
 GGATGCCAGCTGACACTGCAGCAACGTGACTACTGTGCCACTACCTCATCCGGCTGCTGAAGTGAAG
 CGAGACAGCTTCCCAACTTCTGGCCTGCAAGCAGCAGCAGCAGACTGGGACTACTGCGAGCACCTGG
 ATTACGTGAAACGCATGAAGGAGTTTGAACGTGAGCGGCGACTGCTCAAAGGAAGAAGCGAAGGGCACT
 GAAGGAAGCACGGGTAGCCAAGGCCAGGAGAAGGAGAGGTGGGCCCTGAGGTGCCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200842 protein sequence
 Red=Cloning site Green=Tags(s)
 MGAHLTRRYLWDASVEPDPEKIPSPFDLGFPERKERVVMATQQEMMDAQLTLQQRDYCAHYLIRLLKCK
 RDSFPNFLACKHEQHDWDYCEHLDYVKRMKEFERERLLQRKKRRALKEARVAQQGQEGEVGPEVAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


- ACCN:** NM_025843
- ORF Size:** 414 bp
- OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM_025843.3](#), [NP_080119.1](#)
- RefSeq Size:** 533 bp
- RefSeq ORF:** 414 bp
- Locus ID:** 66916


UniProt ID: Q9CR61

Cytogenetics: 8 C2

MW: 16.3 kDa

Gene Summary: Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR200842