

Product datasheet for MR200859

H3f3b (BC092300) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H3f3b (BC092300) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: H3f3b
Synonyms: 9430068D06Rik, H3.3B
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200859 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCCCAACCAAGCAGACCGCTAGGAAGTCCACCGGTGGGAAAGCCCCCGCAAACAGCTGGCCACCA
 AGGCGGCTCGGAAAAGCGCGCCCTCTACCGCGGGGTGAAGAAGCCTCACCGCTACAGGCCAGGGACCGT
 GGCTCTGAGAGAGATCCGTCGTTACCAGAAATCGACTGAGCTGCTCATCCGGAAGCTGCCATTCCAGAGA
 TTGGTGAGGGAGATCGCCAGGATTTCAAACCGACTTGAGGTTTCAAAGTGCAGCCATCGGTGCCCTTC
 AGGAGGCTAGCGAAGCATACCTGGTGGGTTGTTGAAGATACCAATCTGTGTGCCATCCACGCCAAGAG
 AGTCACCATCATGCCAAAGACATCCAGTTGGCTCGCCGGATACGGGGGAGAGAGCTTAAGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200859 protein sequence
 Red=Cloning site Green=Tags(s)
 MARTKQTARKSTGGKAPRKQLATKAARKSAPSTGGVKKPHRYRPGTVALREIRRYQKSTELLIRKLPFQR
 LVREIAQDFKTLRFQSAAGALQEASEAYLVGLFEDTNLCAIHAKRVTIMPKDIQFGSPDTGGESLS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: BC092300

ORF Size: 414 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC092300](#), [AAH92300](#)

RefSeq Size: 1011 bp

RefSeq ORF: 416 bp


Locus ID: 15081

Cytogenetics: 11 E2

MW: 15.3 kDa

Gene Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene contains introns and its mRNA is polyadenylated, unlike most histone genes. The protein encoded by this gene is a replication-independent histone that is a member of the histone H3 family. [provided by RefSeq, Nov 2015]

Product images:



Circular map for MR200859