

## Product datasheet for **MR200859A1V**

### Mouse H3f3b (BC092300) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse H3f3b (BC092300) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** H3f3b  
**Synonyms:** 9430068D06Rik, H3.3B  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200859 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCCGAACCAAGCAGACCGCTAGGAAGTCCACCGGTGGGAAAGCCCCCGCAAACAGCTGGCCACCA  
AGGCGGCTCGGAAAAGCGCGCCCTCTACCGCGGGGTGAAGAAGCCTCACCGCTACAGGCCAGGGACCGT  
GGCTCTGAGAGAGATCCGTCGTTACCAGAAATCGACTGAGCTGCTCATCCGGAAGCTGCCATTCCAGAGA  
TTGGTGAGGGAGATCGCCAGGATTTCAAACCGACTTGAGGTTTCAAAGTGCAGCCATCGGTGCCCTTC  
AGGAGGCTAGCGAAGCATACCTGGTGGGTTGTTTGAAGATACCAATCTGTGTGCCATCCACGCCAAGAG  
AGTCACCATCATGCCAAAGACATCCAGTTTGGCTCGCCGGATACGGGGGAGAGAGCTTAAGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200859 protein sequence  
Red=Cloning site Green=Tags(s)

MARTKQTARKSTGGKAPRKQLATKAARKSAPSTGGVKKPHRYRPGTVALREIRRYQKSTELLIRKLPFQR  
LVREIAQDFKTDLRFQSAALQAEAYLVGLFEDTNLCAIHAKRVTIMPKDIQFGSPDTGGESLS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** BC092300



[View online »](#)

<b>ORF Size:</b>	414 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">BC092300</a></u> , <u><a href="#">AAH92300</a></u>
<b>RefSeq Size:</b>	1011 bp
<b>RefSeq ORF:</b>	416 bp
<b>Locus ID:</b>	15081
<b>Cytogenetics:</b>	11 E2
<b>MW:</b>	15.3 kDa