

Product datasheet for **MR200883A1V**

Mouse Rps13 (BC011192) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Rps13 (BC011192) AAV Particle
Tag: Myc-DDK
Symbol: Rps13
Synonyms: 2700063M04Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200883 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGTGCGATGCACGCTCCCGGAAGGGCCTGTCCCAGTCGGCGCTGCCCTACCGCCGTAGCGTCCCCA
CGTGGCTGAAGTTGACGTCTGACGACGTGAAGGAACAGATTTACAAATTGGCCAAGAAAGGCCTGACTCC
CTCCAGATAGGTGTAATCCTGAGGGACTCACATGGGGTGGCCAGGTCCGTTTTGTGACTGAAATAAA
ATCTTGAGAATCCTCAAGTCAAAGGCCTTGGCCCTGACCTCCCCGAGGATCTCTACCATTTGATTAAGA
AAGCTGTTGCTGTCCGAAAGCACCTTGAGAGGAACAGAAAGGATAAGGATGCTAAATTCGCCTGATTCT
GATAGAGAGCAGAATTCACCGATTGGCTCGATACTATAAGACTAAGCGGGTGTCCCACCCAATTGGAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200883 protein sequence
Red=Cloning site Green=Tags(s)

MGRMHAPGKGLSQSALPYRRSVPTWLKLTSDDVKEQIYKLAKKGLTPSQIGVILRDSHGVAQVRFVTGNK
ILRILKSKGLAPDLPEDLYHLIKKAVAVRKHLELRNRKDKDAKFRLLILIESRIHRLARYYKTKRVLPPNWK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC011192



[View online »](#)

ORF Size:	420 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC011192</u> , <u>AAH11192</u>
RefSeq Size:	1331 bp
RefSeq ORF:	422 bp
Locus ID:	68052
Cytogenetics:	7 F1
MW:	16.1 kDa