

Product datasheet for **MR200887L3V**

Fkbp2 (NM_008020) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fkbp2 (NM_008020) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fkbp2
Synonyms:	13kDa; FKBP-2; mFKBP2; mFKBP13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008020
ORF Size:	423 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200887).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008020.2
RefSeq Size:	661 bp
RefSeq ORF:	423 bp
Locus ID:	14227
UniProt ID:	P45878
Cytogenetics:	19 A
Gene Summary:	PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.[UniProtKB/Swiss-Prot Function]



[View online »](#)