

Product datasheet for **MR200894L4V**

Pfn2 (NM_019410) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Pfn2 (NM_019410) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Pfn2 |
| Synonyms: | Pfn |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_019410 |
| ORF Size: | 423 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR200894). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_019410.2 |
| RefSeq Size: | 2020 bp |
| RefSeq ORF: | 423 bp |
| Locus ID: | 18645 |
| UniProt ID: | Q9JIV2 |
| Cytogenetics: | 3 28.15 cM |
| Gene Summary: | Binds to actin and affects the structure of the cytoskeleton. At high concentrations, profilin prevents the polymerization of actin, whereas it enhances it at low concentrations. By binding to PIP2, it inhibits the formation of IP3 and DG.[UniProtKB/Swiss-Prot Function] |



[View online »](#)