

## Product datasheet for **MR200896L4V**

### 1700028K03Rik (BC106106) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	1700028K03Rik (BC106106) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	1700028K03Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC106106
ORF Size:	420 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200896).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC106106</a>
RefSeq Size:	3105 bp
RefSeq ORF:	422 bp
Locus ID:	76421
Cytogenetics:	5 F



[View online »](#)

**Gene Summary:**

Plays a key role in reinforcing the integrity of the central element of the synaptonemal complex (SC) thereby stabilizing SC, ensuring progression of meiotic prophase I in male and female germ cells (PubMed:30949703). Promotes homologous recombination and crossing-over in meiotic prophase I via its association with SHOC1 (PubMed:30746471). Required for the localization of TEX11 and MSH4 to recombination intermediates (PubMed:30746471). [UniProtKB/Swiss-Prot Function]