

Product datasheet for MR200900

1600014C10Rik (NM_028166) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1600014C10Rik (NM_028166) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 1600014C10Rik
Synonyms: A1428873
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200900 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCCATCATGGTAGATGACATCATGAGGCTTCTGTGCTCCATTTCCAGGAGAGGAAGATGAAGGCGG
 CCGTCAAGCACTCTGAAAGGGTGCCATGGTTGCGGGGGCAATGGCCTTCGTCGGTGGCTTGGTTGGTGG
 CCCACCAGGAATAGCTGTTGGGGGACTGTTGGGGCCTGTTAGGTGCCTGGATGACGAGTGGACAGTTC
 AAGCCCGTTCCTCAGATCTTAATGGAGCTGCCCCCGCTGAACAGCGGAACTTGTAAACGAAGCCATGG
 CCATCATCGGAACTTGGACTGGACAGACGCTGTGCAGCTGACTGCTCTGGTCATGAGCAACCAGGCCAT
 GCAGCAGAGGCTCTTGGCCATGCTAACGACCTATGTCACCAAAGAGCTGCAGGCTGAGATCCGCTATGAG
 GAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200900 protein sequence
 Red=Cloning site Green=Tags(s)

MPIMVDDIMRLLCSISQERKMKA AVKHSKGAMVAGAMAFVGGLVGGPPGIAVGGTVGGLLGAWMTSGQF
 KPVPQILMELPPAEQRKLVNEAMAIIIGNLDWTDVQLTALVMSNQAMQQRLLAML TTYVTKELQAEIRYE
 D

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_028166

ORF Size: 426 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028166.1](#), [NM_028166.2](#), [NM_028166.3](#), [NP_082442.1](#)

RefSeq Size: 3165 bp

RefSeq ORF: 426 bp


Locus ID: 72244

UniProt ID: [Q8WUR0](#)

Cytogenetics: 7 B3

MW: 15.1 kDa

Product images:



Circular map for MR200900