

Product datasheet for **MR200900A1V**

Mouse 1600014C10Rik (NM_028166) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse 1600014C10Rik (NM_028166) AAV Particle
Tag: Myc-DDK
Symbol: 1600014C10Rik
Synonyms: AI428873
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200900 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCCATCATGGTAGATGACATCATGAGGCTTCTGTGCTCCATTTCCAGGAGAGGAAGATGAAGGCGG
 CCGTCAAGCACTCTGAAAGGGTGCCATGGTTGCGGGGCAATGGCCTTCGTCGGTGGCTTGTTGGTGG
 CCCACCAGGAATAGCTGTTGGGGGACTGTTGGGGCCTGTTAGGTGCCTGGATGACGAGTGGACAGTTC
 AAGCCCCTTCAGATCTTAATGGAGCTGCCCCCGCTGAACAGCGGAACTTGTTAACGAAGCCATGG
 CCATCATCGGGAACTTGGACTGGACAGACGCTGTGCAGCTGACTGCTCTGGTCATGAGCAACCAGGCCAT
 GCAGCAGAGGCTCTTGCCATGCTAACGACCTATGTACCAAAGAGCTGCAGGCTGAGATCCGCTATGAG
 GAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200900 protein sequence
 Red=Cloning site Green=Tags(s)

MPIMVDDIMRLLCSISQERKMKAAYKHSGKGAMVAGAMAFVGGLVGGPPGIAVGGTVGGLLGAWMVSGQF
 KPVPQILMELPPAEQRKLVNEAMAIIGNLDWTDVQLTALVMSNQAMQQRLAMLTTYVTKELQAEIRYE
 D

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2



[View online »](#)

ACCN:	NM_028166
ORF Size:	426 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_028166.2</u> , <u>NP_082442.1</u>
RefSeq Size:	3165 bp
RefSeq ORF:	426 bp
Locus ID:	72244
UniProt ID:	<u>Q8WURO</u>
Cytogenetics:	7 B3
MW:	15.1 kDa