

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for MR200907L3V

## C1d (NM\_020558) Mouse Tagged ORF Clone Lentiviral Particle

## **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C1d (NM_020558) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	C1d
Synonyms:	1110036E10Rik; Al875855; SUN-CoR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020558
ORF Size:	426 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200907).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 020558.2, NP 065583.2</u>
RefSeq Size:	2957 bp
RefSeq ORF:	426 bp
Locus ID:	57316
UniProt ID:	<u>O35473</u>
Cytogenetics:	11 A2



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Gene Summary:Plays a role in the recruitment of the RNA exosome complex to pre-rRNA to mediate the 3'-5'<br/>end processing of the 5.8S rRNA; this function may include MPHOSPH6. Can activate PRKDC<br/>not only in the presence of linear DNA but also in the presence of supercoiled DNA. Can<br/>induce apoptosis in a p53/TP53 dependent manner. May regulate the TRAX/TSN complex<br/>formation. Potentiates transcriptional repression by NR1D1 and THRB (By similarity).<br/>[UniProtKB/Swiss-Prot Function]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US