

Product datasheet for MR200908

Plbbp (NM_001039078) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Plbbp (NM_001039078) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Plbbp
Synonyms: 1700024N20Rik; 2200002F22Rik; Prosc
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200908 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTGAGAGGTGGCAGCATGACGGCGGAGCTGGGAGTCGGGTTTCGCGCTGCGGGCCGTGAACGAGCGAG
 TCCAGCAGAGTGTGGCGCGGCCGCGGGACCTCCAGCCATCCAACCCGGCTCGTTGCGGTGAGCAA
 AACAAACCTGCAGACATGGTATCGAGGCCTATGGCCACGGCCAGCGCACTTTGGAGAGAACTATATT
 CTGTCTTCTGTCCCAGATCAAGTGGCACTTCATTGGTCATTTACAGAAACAAAACGTCAACAACTGA
 TGGCCGTCCCAACCTCTCTATGCTGAAACCGTAGACTCGGTGAAGCTGGCAGACAAGGTGAATAGCTC
 CTGGCAGAAGAAAGGTCCTACAGAACCCTGAAGGTGATGGTCCAGATTAAACACCAGCGGAGAGGACAGT
 AAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200908 protein sequence
 Red=Cloning site Green=Tags(s)

MLRGGSMTEALGVGFALRAVNERVQSVARRPRDLPAIQPRLVAVSKTKPADMVEIAYGHGQRTFGENYI
 LSSCPEIKWHFIGHLQKQNVNKLMAVPNLSMLETVDSVKLADKVNSSWQKKGPTTEPLKVMVQINTSGEDS
 K

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online >](#)

Cloning Scheme:



ACCN: NM_001039078

ORF Size: 426 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039078.2](#), [NP_001034167.1](#)

RefSeq Size: 1427 bp

RefSeq ORF: 426 bp


Locus ID: 114863

Cytogenetics: 8 A2

MW: 15.6 kDa

Gene Summary: Pyridoxal 5'-phosphate (PLP)-binding protein, which may be involved in intracellular homeostatic regulation of pyridoxal 5'-phosphate (PLP), the active form of vitamin B6. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR200908