

Product datasheet for **MR200914A1V**

Mouse Med9 (NM_138675) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Med9 (NM_138675) AAV Particle
Tag: Myc-DDK
Symbol: Med9
Synonyms: BC019367; Med25
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200914 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTTCCTCTGGGGTGGCGGGTGGGAGACAGGCCGAAGATACGCTACAGCCACCGCCCGAGCTGCTGC
CGGAGTCCAAGCCGCCACCACCGCTCAGCCTCGCCGGTCGCTGCGCTCCCGCCGCGGGCGCTCCGCG
GCCTCAGTCCCCCGCTGGCGGAAGGAGGAGAATTACTCCTTTTGCCGTTGGTTCACAATGTCATCAA
TGCATGGACAAGGACAGCCCGGATCTCCACCAGGACCTGAACGCCCTCAAACCAAGTCCAGGAGCTGC
GGAAGTCATCGGCACCATGCCCGGCATCCACGTGAGCCCCGAGCAGCAGCAGCAGCTCCACAGCCT
CCGAGAGCAAGTGAGGACCAAGAACGAGCTGCTGCAGAAGTACAAGAGCCTCTGCATGTTTGAGATCCCC
AAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200914 protein sequence
Red=Cloning site Green=Tags(s)

MASSGVAGGRQAEDTLQPPPELLPESKPPPPQPLPVAALPPPAAPRPQSPAGAKEENYSFLPLVHNVIK
CMDKSDPLHQDLNALKTKFQELRKLIGTMPGIHVSPEQQQQLHSLREQVRTKNELLQKYKSLCMFEIP
KD

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2



ACCN:	NM_138675
ORF Size:	426 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_138675.3</u> , <u>NP_619616.2</u>
RefSeq Size:	3013 bp
RefSeq ORF:	429 bp
Locus ID:	192191
UniProt ID:	<u>Q8VCS6</u>
Cytogenetics:	11 B1.3
MW:	15.7 kDa