

Product datasheet for **MR200920**

Fxyd7 (BC061101) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fxyd7 (BC061101) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fxyd7
Synonyms: 1110035I01Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200920 representing BC061101
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCACCCCAACCCAGAGCCCCACAAACGTTCTGAAGAAACAGATCCTTTTTCTATGACTATGCCA
CTGTGCAGACTGTGGGATGACCCTGGCCACTATCATGTTCTGCTGGGATCATCATCCTGAGCAA
GAAGGTGAAGTGCAGGAAGGCGGATCCAGCCCAACATGCAAATCCTGTAAGTCGGAAGTCCCTCCTCA
GCCCTGGAGGTGGCGGTGTGGGATCCTCTGCAGATTTCCACTGCTGTCCCTGGCCAGAGTGCAGGAA
CCAGAGGACCCAGAGAAGGCGGTGGGACCCAGCCTGGCGCCAGGAGTGTGCCCAATCAGCCGCCAG
ACCCACCCCAAGCGCGGAGCCAATGCACCCGCTCTCCCTCCCCAGGCCTTGGCAACAACGATCCTATT
GTACCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200920 representing BC061101
Red=Cloning site Green=Tags(s)

MATPTQSPTNVPEETDPFFYDYATVQTVGMTLATIMFVLGIIIIILSKVKCRKADSSPTCKSCKSELPS
APGGGVWDPLQISTAVPGQSAGTRGPREGGDPAWRQGVCPQSAARPTPRPEPMHPALPPQALATTILF
VP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mm9025_h10.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC061101

ORF Size: 426 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC061101.1](#)

RefSeq Size: 770 bp

RefSeq ORF: 428 bp


Locus ID: 57780

Cytogenetics: 7 B1

MW: 28.2 kDa

Gene Summary: This reference sequence was derived from multiple replicate ESTs and validated by similar human genomic sequence. This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXYD-domain containing ion transport regulator. Transmembrane topology has been established for two family members (FXYP1 and FXYD2), with the N-terminus extracellular and the C-terminus on the cytoplasmic side of the membrane. FXYD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXYD1 (phospholemman), FXYD2 (gamma), FXYD3 (MAT-8), FXYD4 (CHIF), and FXYD5 (RIC) have been shown to induce channel activity in experimental expression systems. This gene product, FXYD7, is novel and has not been characterized as a protein. [RefSeq curation by Kathleen J. Sweadner, Ph.D., sweadner@helix.mgh.harvard.edu., Dec 2000]

Product images:



Circular map for MR200920