

## Product datasheet for **MR200952L3V**

### Pla2g10 (BC028879) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pla2g10 (BC028879) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pla2g10
Synonyms:	PLA2GX, mGXsPLA2, sPLA2-X
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC028879
ORF Size:	432 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200952).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC028879</a> , <a href="#">AAH28879</a>
RefSeq Size:	1011 bp
RefSeq ORF:	434 bp
Locus ID:	26565
Cytogenetics:	16 9.5 cM



[View online »](#)

**Gene Summary:**

This gene encodes a member of the phospholipase A2 family of lipolytic enzymes that hydrolyzes glycerophospholipids to produce free fatty acids and lysophospholipids. The encoded protein undergoes proteolytic processing to generate a calcium-dependent enzyme that plays pivotal roles in the liberation of arachidonic acid from membrane phospholipids leading to the production of various inflammatory lipid mediators, such as prostaglandins. In response to myocardial ischemia/reperfusion, mice lacking the encoded protein display a reduction in myocardial infarct size partly through the suppression of neutrophil cytotoxic activities. Alternative splicing results in multiple transcript variants encoding different isoforms. All of these isoforms may undergo similar processing to generate the mature protein. [provided by RefSeq, Jul 2015]